

Product Datasheet

CD41 antibody (orb434444)

Description	Mouse monoclonal antibody to CD41
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human platelet plasma membranes.
Target	CD41
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.