

Product Datasheet

VON WILLEBRAND FACTOR antibody (orb434448)

Description	Mouse monoclonal antibody to VON WILLEBRAND FACTOR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, Purification
Immunogen	Human Factor VIII complex partially purified from Factor VIII concentrate.
Target	VON WILLEBRAND FACTOR
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.