

Product Datasheet

CYTOCHROME B245 antibody (orb434456)

Description	Mouse monoclonal antibody to CYTOCHROME B245
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr, IP, WB
Immunogen	Cytochrome B245 solubilised in octyl-beta-glucopyranoside.
Target	CYTOCHROME B245
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	LIGHT CHAIN
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.