

Product Datasheet

MHC CLASS II RT1B antibody (Biotin) (orb434475)

Description	Mouse monoclonal antibody against MHC CLASS II RT1B conjugated to Biotin
Species/Host	Mouse
Reactivity	Rat
Conjugation	Biotin
Tested Applications	FC
Immunogen	Rat thymocyte membrane glycoproteins.
Target	MHC CLASS II RT1B
Preservatives	0.09% Sodium Azide1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.