

Product Datasheet

MHC CLASS II RT1B antibody (orb434479)

Description	Mouse monoclonal antibody to MHC CLASS II RT1B
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IHC-P
Immunogen	Rat thymocyte membrane glycoproteins.
Target	MHC CLASS II RT1B
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.