


Product Datasheet

ADCY10 antibody (orb4345)

Description	ADCY10 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH conjugated synthetic peptide derived from human ADCY10 (1101-1200/1610 aa)
Target	ADCY10
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	187 kDa
Uniprot ID	Q96PN6
Dilution Range	WB=1:500-2000, IF=1:200, IHC-P=1:100-500, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...