

Product Datasheet

P2X7 antibody (orb434500)

Description	Rat monoclonal antibody to P2X7
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF
Immunogen	A P2X7-expression construct and a final boost with P2X7-transfected HEK cells.
Target	P2X7
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.