

Product Datasheet

GFAP antibody (orb434511)

Description	Mouse monoclonal antibody to GFAP
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	IH-Fr, IHC-P, RIA, WB
Immunogen	Bovine spinal cord homogenate.
Target	GFAP
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.