

Product Datasheet

CLOSTRIDIUM DIFFICILE TOXIN B antibody (orb434515)

Description	Mouse monoclonal antibody to CLOSTRIDIUM DIFFICILE TOXIN B
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Target	CLOSTRIDIUM DIFFICILE TOXIN B
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.