

Product Datasheet

GAPDH antibody (HRP) (orb434518)

Description	Mouse monoclonal antibody against GAPDH conjugated to HRP
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, WB
Immunogen	Rabbit muscle GAPDH.
Target	GAPDH
Preservatives	0.01% Thiomersal. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.