

Product Datasheet

RHODOPSIN antibody (orb434521)

Description	Mouse monoclonal antibody to RHODOPSIN
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Immunogen	Bovine rhodopsin.
Target	RHODOPSIN
Preservatives	0.09% Sodium Azide (NaN ₃)0.01% Bovine Serum Albumin50% Glycerol. 10mM HEPES pH 7.5, 150mM NaCl
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.