

Product Datasheet

NITROTYROSINE antibody (orb434539)

Description	Mouse monoclonal antibody to NITROTYROSINE
Species/Host	Mouse
Reactivity	Chemical
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Protein-bound Nitrotyrosine.
Target	NITROTYROSINE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.