

Product Datasheet

DYKDDDDK TAG antibody (orb434541)

Description	Rat monoclonal antibody to DYKDDDDK TAG
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF, IHC-P, IP, WB
Target	DYKDDDDK TAG
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.