

Product Datasheet

HLA DP DQ DR antibody (orb434545)

Description	Mouse monoclonal antibody to HLA DP DQ DR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FACS, IHC-P
Immunogen	Human HLA Class II (DP, DQ, DR).
Target	HLA DP DQ DR
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.