

## Product Datasheet

### F4/80 antibody (orb434589)

|                            |                                                                     |
|----------------------------|---------------------------------------------------------------------|
| <b>Description</b>         | Rat monoclonal antibody to F4/80                                    |
| <b>Species/Host</b>        | Rat                                                                 |
| <b>Reactivity</b>          | Mouse                                                               |
| <b>Conjugation</b>         | Unconjugated                                                        |
| <b>Tested Applications</b> | FACS, ICC, IF, IHC-Fr, IHC-P, IP                                    |
| <b>Immunogen</b>           | Thioglycollate stimulated peritoneal macrophages from C57BL/6 mice. |
| <b>Target</b>              | F4/80                                                               |
| <b>Preservatives</b>       | 0.09% Sodium Azide. None present                                    |
| <b>Form/Appearance</b>     | Tissue Culture Supernatant - liquid                                 |
| <b>Storage</b>             | +4°C, -20°C if preferred                                            |
| <b>Note</b>                | For research use only                                               |
| <b>Isotype</b>             | IgG2b                                                               |
| <b>Clonality</b>           | Monoclonal                                                          |
| <b>Purity</b>              | S/N                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                     |