

Product Datasheet

F4/80 antibody (orb434604)

Description	Rat monoclonal antibody to F4/80
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IHC-P, IP, WB
Immunogen	Thioglycollate stimulated peritoneal macrophages from C57BL/6 mice.
Target	F4/80
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.