

Product Datasheet

CD18 antibody (RPE) (orb434619)

Description	Rat monoclonal antibody against CD18 conjugated to RPE
Species/Host	Rat
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Human neutrophils
Target	CD18
Preservatives	0.09% Sodium Azide1%Bovine Serum Albumin5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.