

Product Datasheet

CD235a antibody (FITC) (orb434622)

Description	Rat monoclonal antibody against CD235a conjugated to FITC
Species/Host	Rat
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Target	CD235a
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Ig Fraction conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.