

Product Datasheet

CD235a antibody (orb434623)

Description	Rat monoclonal antibody to CD235a
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Cytotoxic Assay, FC, IH-Fr, IHC-P
Target	CD235a
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.