

Product Datasheet

MACROPHAGES/MONOCYTES antibody (orb434635)

Description	Rat monoclonal antibody to MACROPHAGES/MONOCYTES
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr
Immunogen	Mouse lymph node stroma.
Target	MACROPHAGES/MONOCYTES
Preservatives	0.1% Sodium Azide (NaN ₃).. Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.