

Product Datasheet

MHC CLASS I RT1A antibody (orb434641)

Description	Mouse monoclonal antibody to MHC CLASS I RT1A
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr, IP
Immunogen	Rat spleen cell glycoproteins
Target	MHC CLASS I RT1A
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.