

Product Datasheet

CD45RC antibody (orb434650)

Description	Mouse monoclonal antibody to CD45RC
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P
Immunogen	Phytohaemagglutinin (PHA) -activated rat lymphocytes
Target	CD45RC
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.