

Product Datasheet

CD45RC antibody (orb434652)

Description	Mouse monoclonal antibody to CD45RC
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P
Immunogen	Phytohaemagglutinin (PHA) -activated rat lymphocytes
Target	CD45RC
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.