

## Product Datasheet

### CD45RC antibody (orb434653)

|                            |                                                      |
|----------------------------|------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CD45RC                  |
| <b>Species/Host</b>        | Mouse                                                |
| <b>Reactivity</b>          | Rat                                                  |
| <b>Conjugation</b>         | Unconjugated                                         |
| <b>Tested Applications</b> | FC, IH-Fr, IHC-P                                     |
| <b>Immunogen</b>           | Phytohaemagglutinin (PHA) -activated rat lymphocytes |
| <b>Target</b>              | CD45RC                                               |
| <b>Preservatives</b>       | 0.09% Sodium Azide. Phosphate buffered saline        |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                |
| <b>Concentration</b>       | 1.0mg/ml                                             |
| <b>Storage</b>             | +4°C, -20°C if preferred                             |
| <b>Note</b>                | For research use only                                |
| <b>Isotype</b>             | IgG1                                                 |
| <b>Clonality</b>           | Monoclonal                                           |
| <b>Purity</b>              | Purified                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                      |