

Product Datasheet

CD35 antibody (orb434681)

Description	Mouse monoclonal antibody to CD35
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P
Immunogen	Human acute monocytic leukaemia cells.
Target	CD35
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.