

Product Datasheet

APOLIPOPROTEIN E antibody (Biotin) (orb434697)

Description	Mouse monoclonal antibody against APOLIPOPROTEIN E conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Purified ApoHDL fraction.
Target	APOLIPOPROTEIN E
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.