

Product Datasheet

LIPOPROTEIN LIPASE antibody (orb434701)

Description	Mouse monoclonal antibody to LIPOPROTEIN LIPASE
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Purified bovine milk lipoprotein lipase.
Target	LIPOPROTEIN LIPASE
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.