

Product Datasheet

LIPOPROTEIN LIPASE antibody (Biotin) (orb434702)

Description	Mouse monoclonal antibody against LIPOPROTEIN LIPASE conjugated to Biotin
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA, IP, WB
Immunogen	Purified bovine milk lipoprotein lipase.
Target	LIPOPROTEIN LIPASE
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.