

Product Datasheet

LIPOPROTEIN LIPASE antibody (orb434703)

Description	Mouse monoclonal antibody to LIPOPROTEIN LIPASE
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA, FA, IP, WB
Immunogen	Purified bovine milk lipoprotein lipase.
Target	LIPOPROTEIN LIPASE
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.