

Product Datasheet

LIPOPROTEIN LIPASE antibody (HRP) (orb434705)

Description	Mouse monoclonal antibody against LIPOPROTEIN LIPASE conjugated to HRP
Species/Host	Mouse
Reactivity	Bovine
Conjugation	HRP
Tested Applications	ELISA, IP, WB
Immunogen	Purified bovine milk lipoprotein lipase.
Target	LIPOPROTEIN LIPASE
Preservatives	0.01% Thiomersal. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.