

Product Datasheet

CHROMOGRANIN A antibody (orb434749)

Description	Mouse monoclonal antibody to CHROMOGRANIN A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, IHC-P, WB
Immunogen	Human pheochromocytoma tissue.
Target	CHROMOGRANIN A
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.