

Product Datasheet

NOTCH 2 antibody (orb434757)

Description	Hamster monoclonal antibody to NOTCH 2
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Mouse Notch 2-Fc fusion protein.
Target	NOTCH 2
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.