

Product Datasheet

DELTA-LIKE PROTEIN 1 antibody (RPE) (orb434762)

Description	Hamster monoclonal antibody against DELTA-LIKE PROTEIN 1 conjugated to RPE
Species/Host	Hamster
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	DLL1-expressing CHO cells.
Target	DELTA-LIKE PROTEIN 1
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.