

Product Datasheet

DELTA-LIKE PROTEIN 4 antibody (orb434764)

Description	Hamster monoclonal antibody to DELTA-LIKE PROTEIN 4
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FACS, IHC-Fr, IHC-P
Immunogen	Recombinant mouse DLL4.
Target	DELTA-LIKE PROTEIN 4
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.