

Product Datasheet

JAGGED2 antibody (orb434768)

Description	Hamster monoclonal antibody to JAGGED2
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P
Immunogen	Jagged2-expressing CHO cells.
Target	JAGGED2
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.