



Product Datasheet

CD68 antibody (RPE) (orb434771)


Description	Mouse monoclonal antibody against CD68 conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Subcellular fraction of human alveolar macrophages
Target	CD68
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 5% Sucrose.
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	Prior to reconstitution store at +4°C. After reconstitution store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend micro-centrifugation before use.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
MW	110 kDa
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human to...



Immunohistochemical staining of human br...



Flow cytometric analysis of Human CD14 u...