

Product Datasheet

VEGF antibody (orb434807)

Description	Mouse monoclonal antibody to VEGF
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, WB
Target	VEGF
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.