

Product Datasheet

EOSINOPHIL MAJOR BASIC PROTEIN antibody (orb434808)

Description	Mouse monoclonal antibody to EOSINOPHIL MAJOR BASIC PROTEIN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, IHC-P
Target	EOSINOPHIL MAJOR BASIC PROTEIN
Preservatives	0.02% Sodium Azide (NaN ₃)0.1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.