

Product Datasheet

MHC CLASS I antibody (orb434810)

Description	Mouse monoclonal antibody to MHC CLASS I
Species/Host	Mouse
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP, WB
Target	MHC CLASS I
Preservatives	0.1% Sodium Azide (NaN ₃). Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.