

Product Datasheet

HIGH MOBILITY GROUP PROTEIN B1 antibody (orb434814)

Description	Hamster monoclonal antibody to HIGH MOBILITY GROUP PROTEIN B1
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, WB
Immunogen	Recombinant mouse HMGB1.
Target	HIGH MOBILITY GROUP PROTEIN B1
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.