

Product Datasheet

TCR ALPHA/BETA antibody (orb434818)

Description	Mouse monoclonal antibody to TCR ALPHA/BETA
Species/Host	Mouse
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr, IP
Target	TCR ALPHA/BETA
Preservatives	0.1% Sodium Azide (NaN ₃). Borate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	Vb1
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.