

Product Datasheet

TCR GAMMA/DELTA antibody (orb434820)

Description	Mouse monoclonal antibody to TCR GAMMA/DELTA
Species/Host	Mouse
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Target	TCR GAMMA/DELTA
Preservatives	0.1% Sodium Azide (NaN ₃). Borate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.