

Product Datasheet

CD28 antibody (orb434825)

Description	Mouse monoclonal antibody to CD28
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	FA, FACS
Immunogen	Plasmid DNA encoding the bovine CD28 gene.
Target	CD28
Preservatives	None present. Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.