

Product Datasheet

TIM-1 antibody (FITC) (orb434846)

Description	Rat monoclonal antibody against TIM-1 conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Full-length mouse Tim-1-Ig, containing both the IgV and mucin domains of Tim-1.
Target	TIM-1
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.