

Product Datasheet

TIM-1 antibody (orb434847)

Description	Rat monoclonal antibody to TIM-1
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF
Immunogen	Full-length mouse Tim-1-Ig, containing both the IgV and mucin domains of Tim-1.
Target	TIM-1
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.