

Product Datasheet

TIM-1 antibody (RPE) (orb434848)

Description	Rat monoclonal antibody against TIM-1 conjugated to RPE
Species/Host	Rat
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	Full-length mouse Tim-1-Ig, containing both the IgV and mucin domains of Tim-1.
Target	TIM-1
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.