

Product Datasheet

TIM-2 antibody (orb434851)

Description	Rat monoclonal antibody to TIM-2
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Tim-2-Ig, consisting of amino acids 1-230 of mouse Tim-2, coupled with the Fc portion of mouse IgG2a.
Target	TIM-2
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.