

Product Datasheet

TIM-3 antibody (FITC) (orb434854)

Description	Rat monoclonal antibody against TIM-3 conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Tim-3-Ig, consisting of amino acids 1-191 of mouse Tim-3, coupled with the Fc portion of mouse IgG2a.
Target	TIM-3
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.