

## Product Datasheet

### CD283 antibody (orb434874)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | Rat monoclonal antibody to CD283                   |
| <b>Species/Host</b>        | Rat                                                |
| <b>Reactivity</b>          | Mouse                                              |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | FACS, ICC, IF                                      |
| <b>Immunogen</b>           | Purified murine TLR3 extracellular domain protein. |
| <b>Target</b>              | CD283                                              |
| <b>Preservatives</b>       | 0.09% Sodium Azide.. Phosphate buffered saline     |
| <b>Form/Appearance</b>     | Purified IgG - liquid                              |
| <b>Concentration</b>       | 1.0mg/ml                                           |
| <b>Storage</b>             | +4°C, -20°C if preferred                           |
| <b>Note</b>                | For research use only                              |
| <b>Isotype</b>             | IgG2a                                              |
| <b>Clonality</b>           | Monoclonal                                         |
| <b>Purity</b>              | Purified                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |