

Product Datasheet

CD32 antibody (RPE) (orb434882)

Description	Mouse monoclonal antibody against CD32 conjugated to RPE
Species/Host	Mouse
Reactivity	Bovine
Conjugation	RPE
Tested Applications	FC
Immunogen	Bovine FcγRII-transfected COS7 cells.
Target	CD32
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.